

For Research Use Only

# SPP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60558-6-PBS



## Basic Information

Catalog Number:

60558-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16823

GenBank Accession Number:

BC008938

GeneID (NCBI):

81502

UNIPROT ID:

Q8TCT9

Full Name:

histocompatibility (minor) 13

Calculated MW:

377 aa, 41 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H9C9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50787-4, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-5-PBS. Detection antibody: 60558-6-PBS. Standard:Ag16823. Range: 0.391-100 ng/mL.