

For Research Use Only

PTRF Monoclonal antibody

Catalog Number:60562-4-Ig



Basic Information

Catalog Number:

60562-4-Ig

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG13495

GenBank Accession Number:

BC008849

GeneID (NCBI):

284119

UNIPROT ID:

Q6NZI2

Full Name:

polymerase I and transcript release factor

Calculated MW:

43 kDa

Observed MW:

50-55 kDa

Purification Method:

Protein A purification

CloneNo.:

3B9F11

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat, pig, rabbit

Positive Controls:

WB : HaCaT cells, A431 cells, A549 cells, HeLa cells, K-562 cells, rabbit colon tissue, rat colon tissue, mouse colon tissue, pig adipose tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60562-4-Ig (PTRF antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.