

For Research Use Only

PTRF Monoclonal antibody, PBS Only

Catalog Number:60562-4-PBS



Basic Information

Catalog Number:

60562-4-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG13495

GenBank Accession Number:

BC008849

GeneID (NCBI):

284119

UNIPROT ID:

Q6NZI2

Full Name:

polymerase I and transcript release factor

Calculated MW:

43 kDa

Observed MW:

50-55 kDa

Purification Method:

Protein A purification

CloneNo.:

3B9F11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat, pig, rabbit

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60562-4-Ig (PTRF antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60562-4-PBS in a different storage buffer formulation.