

For Research Use Only

Integrin alpha-7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60564-2-PBS



Basic Information

Catalog Number:

60564-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19777

GenBank Accession Number:

BC050280

GeneID (NCBI):

3679

UNIPROT ID:

Q13683

Full Name:

integrin, alpha 7

Calculated MW:

1181 aa, 129 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D11H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50795-1, Integrin alpha-7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60564-1-PBS. Detection antibody: 60564-2-PBS. Standard:Ag19777. Range: 0.098-6.25 ng/mL



Cytometric bead array standard curve of MP50795-2, Integrin alpha-7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60564-3-PBS. Detection antibody: 60564-2-PBS. Standard:Ag19777. Range: 1.563-50 ng/mL