

For Research Use Only

SCD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60569-1-PBS



Basic Information

Catalog Number:

60569-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29167

GenBank Accession Number:

BC005807

GeneID (NCBI):

6319

UNIPROT ID:

O00767

Full Name:

stearoyl-CoA desaturase (delta-9-desaturase)

Calculated MW:

355 aa, 41 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F9B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50806-1, SCD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60569-1-PBS. Detection antibody: 60569-2-PBS. Standard:Ag29167. Range: 6.25-100 ng/mL