

For Research Use Only

CYP2C9 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60576-5-PBS



Basic Information

Catalog Number:

60576-5-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9566

GenBank Accession Number:

BC020754

GeneID (NCBI):

1559

UNIPROT ID:

P11712

Full Name:

cytochrome P450, family 2, subfamily
C, polypeptide 9

Calculated MW:

162aa, 18 kDa; 490aa, 56 kDa

Purification Method:

Protein A purification

CloneNo.:

2B1G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50815-4, CYP2C9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60576-5-PBS. Detection antibody: 60576-2-PBS. Standard: Ag9566. Range: 1.563-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 60576-5-Ig against Human CYP2C9 was performed. The affinity constant is below 1 pM.