

For Research Use Only

RSAD2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60591-2-PBS



Basic Information

Catalog Number:

60591-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG15174

GenBank Accession Number:

BC017969

GeneID (NCBI):

91543

UNIPROT ID:

Q8WXG1

Full Name:

radical S-adenosyl methionine
domain containing 2

Calculated MW:

361 aa, 42 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G6G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50840-1, RSAD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60591-1-PBS. Detection antibody: 60591-2-PBS. Standard:Ag15174. Range: 6.25-100 ng/mL