

For Research Use Only

SERCA2/ATP2A2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60593-4-PBS



Basic Information

Catalog Number:

60593-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26307

GenBank Accession Number:

BC035588

GeneID (NCBI):

488

UNIPROT ID:

P16615

Full Name:

ATPase, Ca⁺⁺ transporting, cardiac
muscle, slow twitch 2

Calculated MW:

1042 aa, 115 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F10C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50843-2, SERCA2,ATP2A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60593-3-PBS. Detection antibody: 60593-4-PBS. Standard:Ag5052. Range: 1.563-100 ng/mL