

For Research Use Only

# MARCH5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60603-2-PBS



## Basic Information

Catalog Number:

60603-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8286

GenBank Accession Number:

BC015480

GeneID (NCBI):

54708

UNIPROT ID:

Q9NX47

Full Name:

membrane-associated ring finger  
(C3HC4) 5

Calculated MW:

278 aa, 31 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4G3A5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50859-1, MARCH5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60603-1-PBS. Detection antibody: 60603-2-PBS. Standard:Ag8286. Range: 0.098-25 ng/mL