

For Research Use Only

PSEN2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60609-1-PBS



Basic Information

Catalog Number:

60609-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33489

GenBank Accession Number:

BC006365

GeneID (NCBI):

5664

UNIPROT ID:

P49810

Full Name:

presenilin 2 (Alzheimer disease 4)

Calculated MW:

415 aa, 46 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3F1A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50866-1, PSEN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60609-1-PBS. Detection antibody: 60609-2-PBS. Standard:Ag33489. Range: 3.125-100 ng/mL.