

For Research Use Only

CDO1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60610-1-PBS



Basic Information

Catalog Number:

60610-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3294

GenBank Accession Number:

BC024241

GeneID (NCBI):

1036

UNIPROT ID:

Q16878

Full Name:

cysteine dioxygenase, type I

Calculated MW:

200 aa, 23 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D2E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array standard curve of MP50867-1, CDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60610-1-PBS. Detection antibody: 60610-2-PBS. Standard:Ag3294. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50867-2, CDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60610-1-PBS. Detection antibody: 60610-3-PBS. Standard:Ag3294. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50867-3, CDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60610-1-PBS. Detection antibody: 60610-4-PBS. Standard:Ag3294. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50867-4, CDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60610-5-PBS. Detection antibody: 60610-1-PBS. Standard:Ag3294. Range: 0.098-100 ng/mL