

For Research Use Only

RDH5 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60617-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60617-1-PBS	BC028298	Protein G Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	5959	1H3H12
Source:	UNIPROT ID:	
Mouse	Q92781	
Isotype:	Full Name:	
IgG1	retinol dehydrogenase 5 (11-cis/9-cis)	
Immunogen Catalog Number:	Calculated MW:	
AG30513	318 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50879-1, RDH5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60617-1-PBS. Detection antibody: 60617-2-PBS. Standard:Ag30513. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50879-2, RDH5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60617-3-PBS. Detection antibody: 60617-1-PBS. Standard:Ag30513. Range: 0.098-100 ng/mL