

For Research Use Only

# RDH5 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60617-2-PBS



## Basic Information

Catalog Number:

60617-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG30513

GenBank Accession Number:

BC028298

GeneID (NCBI):

5959

UNIPROT ID:

Q92781

Full Name:

retinol dehydrogenase 5 (11-cis/9-cis)

Calculated MW:

318 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

2D5E4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50879-1, RDH5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60617-1-PBS. Detection antibody: 60617-2-PBS. Standard:Ag30513. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50879-3, RDH5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60617-2-PBS. Detection antibody: 60617-4-PBS. Standard:Ag30513. Range: 0.098-100 ng/mL