

For Research Use Only

CYP19A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60624-3-PBS



Basic Information

Catalog Number:

60624-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9896

GenBank Accession Number:

BC022896

GeneID (NCBI):

1588

UNIPROT ID:

P11511

Full Name:

cytochrome P450, family 19, subfamily A, polypeptide 1

Calculated MW:

20 kDa, 25 kDa, 58 kDa

Purification Method:

Protein G purification

CloneNo.:

1A6G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50890-2, CYP19A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60624-3-PBS. Detection antibody: 60624-4-PBS. Standard:Ag9896. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50890-3, CYP19A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60624-3-PBS. Detection antibody: 60624-5-PBS. Standard:Ag9896. Range: 0.098-100 ng/mL.