

For Research Use Only

CLEC2D Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60627-3-PBS



Basic Information

Catalog Number:

60627-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG31670

GenBank Accession Number:

BC019883

GeneID (NCBI):

29121

UNIPROT ID:

Q9UHP7

Full Name:

C-type lectin domain family 2,
member D

Calculated MW:

154 aa, 18 kDa

Purification Method:

Protein A purification

CloneNo.:

1C11B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50893-2, CLEC2D Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60627-3-PBS. Detection antibody: 60627-2-PBS. Standard:Ag31670. Range: 0.781-100 ng/mL.