

For Research Use Only

# RHBDF2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60630-3-PBS



## Basic Information

Catalog Number:

60630-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19855

GenBank Accession Number:

BC016034

GeneID (NCBI):

79651

UNIPROT ID:

Q6PJF5

Full Name:

rhomboid 5 homolog 2 (Drosophila)

Calculated MW:

856 aa, 97 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A12F11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50897-1, RHBDP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60630-2-PBS. Detection antibody: 60630-3-PBS. Standard:Ag19855. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50897-3, RHBDP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60630-5-PBS. Detection antibody: 60630-3-PBS. Standard:Ag19855. Range: 0.098-100 ng/mL