

For Research Use Only

MZB1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60633-3-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60633-3-PBS	BC021275	Protein G Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	51237	1G10C7
Source:	UNIPROT ID:	
Mouse	Q8WU39	
Isotype:	Full Name:	
IgG1	hypothetical protein MGC29506	
Immunogen Catalog Number:	Calculated MW:	
AG17912	189 aa, 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50905-2, MZB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60633-1-PBS. Detection antibody: 60633-3-PBS. Standard:Ag17912. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50905-3, MZB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60633-2-PBS. Detection antibody: 60633-3-PBS. Standard:Ag17912. Range: 0.098-100 ng/mL