

For Research Use Only

# ABCG1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60637-1-PBS



## Basic Information

Catalog Number:

60637-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17329

GenBank Accession Number:

BC029158

GeneID (NCBI):

9619

UNIPROT ID:

P45844

Full Name:

ATP-binding cassette, sub-family G (WHITE), member 1

Calculated MW:

678 aa, 76 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E4G2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50911-1, ABCG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60637-1-PBS. Detection antibody: 60637-2-PBS. Standard:Ag17329. Range: 0.391-100 ng/mL.