

For Research Use Only

SPAG5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60644-3-PBS



Basic Information

Catalog Number:

60644-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6550

GenBank Accession Number:

BC000322

GeneID (NCBI):

10615

UNIPROT ID:

Q96R06

Full Name:

sperm associated antigen 5

Calculated MW:

134 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1D10D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50922-2, SPAG5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60644-3-PBS. Detection antibody: 60644-4-PBS. Standard:Ag6550. Range: 3.125-200 ng/mL.