

For Research Use Only

AGPAT3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60649-1-PBS



Basic Information

Catalog Number:

60649-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG22855

GenBank Accession Number:

BC011971

GeneID (NCBI):

56894

UNIPROT ID:

Q9NRZ7

Full Name:

1-acylglycerol-3-phosphate O-acyltransferase 3

Calculated MW:

376 aa, 43 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1E3D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50928-1, AGPAT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60649-1-PBS. Detection antibody: 60649-2-PBS. Standard:Ag22855. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50928-2, AGPAT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60649-1-PBS. Detection antibody: 60649-3-PBS. Standard:Ag22855. Range: 0.195-100 ng/mL