

For Research Use Only

Phospholamban Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60659-2-PBS



Basic Information

Catalog Number:

60659-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG8632

GenBank Accession Number:

BC005269

GeneID (NCBI):

5350

UNIPROT ID:

P26678

Full Name:

phospholamban

Calculated MW:

52 aa, 6 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

4C9E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50943-1, PLN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60659-1-PBS. Detection antibody: 60659-2-PBS. Standard:Ag8632. Range: 0.098-100 ng/mL.