

For Research Use Only

SEPT12 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60661-4-PBS



Basic Information

Catalog Number:

60661-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4644

GenBank Accession Number:

BC035619

GeneID (NCBI):

124404

UNIPROT ID:

Q8IYM1

Full Name:

septin 12

Calculated MW:

358 aa, 41 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D9B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50944-2, SEPT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60661-3-PBS. Detection antibody: 60661-4-PBS. Standard:Ag4644. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50944-3, SEPT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60661-5-PBS. Detection antibody: 60661-4-PBS. Standard:Ag4644. Range: 0.391-100 ng/mL