

For Research Use Only

# SLC43A1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60667-4-PBS



## Basic Information

Catalog Number:

60667-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31817

GenBank Accession Number:

BC001639

GeneID (NCBI):

8501

UNIPROT ID:

O75387

Full Name:

solute carrier family 43, member 1

Calculated MW:

61 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C6D4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50952-2, SLC43A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60667-3-PBS. Detection antibody: 60667-4-PBS. Standard:Ag31817. Range: 1.563-200 ng/mL.