

For Research Use Only

QSOX1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60676-1-PBS



Basic Information

Catalog Number: 60676-1-PBS	GenBank Accession Number: BC017692	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 5768	CloneNo.: 1D12G10
Source: Mouse	UNIPROT ID: O00391	
Isotype: IgG1	Full Name: quiescin Q6 sulfhydryl oxidase 1	
Immunogen Catalog Number: AG17036	Calculated MW: 747 aa, 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP50971-1, Human QSOX1/Quiescin Q6 Monoclonal Matched Antibody Pair - PBS only. 60676-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17036. 60676-2-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



Serum of thirteen individual healthy human donors was measured. The QSOX1/Quiescin Q6 concentration of detected samples was determined to be 26.07 ng/mL with a range of 7.70-126.24 ng/mL



Cytometric bead array standard curve of MP50971-1, QSOX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60676-1-PBS. Detection antibody: 60676-2-PBS. Standard:Ag17036. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50971-2, QSOX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60676-1-PBS. Detection antibody: 60676-3-PBS. Standard:Ag17036. Range: 0.391-100 ng/mL