

For Research Use Only

# STIM2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60684-3-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60684-3-PBS	BC146661	Protein G Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	57620	2C10H11
Source:	UNIPROT ID:	
Mouse	Q9P246	
Isotype:	Full Name:	
IgG1	stromal interaction molecule 2	
Immunogen Catalog Number:	Calculated MW:	
AG15530	754 aa, 85 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50986-2, STIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60684-1-PBS. Detection antibody: 60684-3-PBS. Standard:Ag15530. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50986-3, STIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60684-2-PBS. Detection antibody: 60684-3-PBS. Standard:Ag15530. Range: 0.098-100 ng/mL