

For Research Use Only

EXT2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60686-2-PBS



Basic Information

Catalog Number: 60686-2-PBS	GenBank Accession Number: BC010058	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2132	CloneNo.: 1D1E11
Source: Mouse	UNIPROT ID: Q93063	
Isotype: IgG1	Full Name: exostoses (multiple) 2	
Immunogen Catalog Number: AG1748	Calculated MW: 718 aa, 82 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50989-1, EXT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60686-1-PBS. Detection antibody: 60686-2-PBS. Standard:Ag1748. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50989-3, EXT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60686-4-PBS. Detection antibody: 60686-2-PBS. Standard:Ag1748. Range: 3.125-100 ng/mL