

For Research Use Only

BBS7 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60687-1-PBS



Basic Information

Catalog Number:

60687-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4601

GenBank Accession Number:

BC032691

GeneID (NCBI):

55212

UNIPROT ID:

Q8IWZ6

Full Name:

Bardet-Biedl syndrome 7

Calculated MW:

715 aa, 80 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C5G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50991-1, BBS7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60687-1-PBS. Detection antibody: 60687-2-PBS. Standard:Ag4601. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50991-2, BBS7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60687-1-PBS. Detection antibody: 60687-3-PBS. Standard:Ag4601. Range: 0.098-12.5 ng/mL.