

For Research Use Only

# GSDMD Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60693-2-PBS



## Basic Information

Catalog Number:

60693-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33638

GenBank Accession Number:

BC008904

GeneID (NCBI):

79792

UNIPROT ID:

P57764

Full Name:

gasdermin D

Calculated MW:

484 aa, 53 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F12C3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51002-1, GSDMD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60693-1-PBS. Detection antibody: 60693-2-PBS. Standard:Ag16827. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51002-3, GSDMD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60693-5-PBS. Detection antibody: 60693-2-PBS. Standard:Ag16827. Range: 1.563-100 ng/mL