

For Research Use Only

KDEL^R1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60694-2-PBS



Basic Information

Catalog Number:

60694-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG26298

GenBank Accession Number:

BC018778

GeneID (NCBI):

10945

UNIPROT ID:

P24390

Full Name:

KDEL (Lys-Asp-Glu-Leu) endoplasmic
reticulum protein retention receptor 1

Calculated MW:

212 aa, 25 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3F10E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51003-1, KDELR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60694-1-PBS. Detection antibody: 60694-2-PBS. Standard:Ag26298. Range: 0.391-100 ng/mL.