

For Research Use Only

NR4A3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60697-1-PBS



Basic Information

Catalog Number:

60697-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28956

GenBank Accession Number:

NM_006981

GeneID (NCBI):

8013

UNIPROT ID:

Q92570

Full Name:

nuclear receptor subfamily 4, group A, member 3

Calculated MW:

68 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G4G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51007-1, NR4A3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60697-1-PBS. Detection antibody: 60697-2-PBS. Standard:Ag28956. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP51007-2, NR4A3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60697-1-PBS. Detection antibody: 60697-3-PBS. Standard:Ag28956. Range: 0.098-100 ng/mL.