

For Research Use Only

# LTBP2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60701-5-PBS



## Basic Information

Catalog Number:

60701-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12718

GenBank Accession Number:

BC078659

GeneID (NCBI):

4053

UNIPROT ID:

Q14767

Full Name:

latent transforming growth factor beta  
binding protein 2

Calculated MW:

1821 aa, 195 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4H5H12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51013-3, LTBP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60701-5-PBS. Detection antibody: 60701-2-PBS. Standard:Ag12718. Range: 0.391-25 ng/mL