

For Research Use Only

DST Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60703-1-PBS



Basic Information

Catalog Number:

60703-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28858

GenBank Accession Number:

NM_001144769

GeneID (NCBI):

667

UNIPROT ID:

Q03001

Full Name:

dystonin

Calculated MW:

861 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E5D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP51015-1, DST Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60703-1-PBS. Detection antibody: 60703-2-PBS. Standard:Ag28858. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51015-2, DST Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60703-1-PBS. Detection antibody: 60703-3-PBS. Standard:Ag28858. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51015-3, DST Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60703-1-PBS. Detection antibody: 60703-4-PBS. Standard:Ag28858. Range: 0.098-100 ng/mL