

For Research Use Only

ATF4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60706-2-PBS



Basic Information

Catalog Number:

60706-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16665

GenBank Accession Number:

BC022088

GeneID (NCBI):

468

UNIPROT ID:

P18848

Full Name:

activating transcription factor 4 (tax-responsive enhancer element B67)

Calculated MW:

39 kDa

Purification Method:

Protein G purification

CloneNo.:

3G3A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51019-1, ATF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60706-1-PBS. Detection antibody: 60706-2-PBS. Standard:Ag16665. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51019-4, ATF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60035-3-PBS. Detection antibody: 60706-2-PBS. Standard:Ag16665. Range: 1.563-100 ng/mL