

For Research Use Only

NELL1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60707-1-PBS



Basic Information

Catalog Number:

60707-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27773

GenBank Accession Number:

BC069674

GeneID (NCBI):

4745

UNIPROT ID:

Q92832

Full Name:

NEL-like 1 (chicken)

Calculated MW:

810 aa, 90 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A9H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51021-1, NEU11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60707-1-PBS. Detection antibody: 60707-2-PBS. Standard:Ag27773. Range: 0.781-200 ng/mL.