

For Research Use Only

# IFT20 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60714-1-PBS



## Basic Information

Catalog Number:

60714-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4718

GenBank Accession Number:

BC038094

GeneID (NCBI):

90410

UNIPROT ID:

Q8IY31

Full Name:

intraflagellar transport 20 homolog  
(Chlamydomonas)

Calculated MW:

15 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5D12C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array sample test of MP51743-1, IFT20 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60765-2-PBS. Detection antibody: 60714-1-PBS.



Cytometric bead array standard curve of MP51743-1, IFT20 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60765-2-PBS. Detection antibody: 60714-1-PBS. Standard:Ag4718. Range: 0.781-100 ng/mL.