

For Research Use Only

TNFRSF11B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60715-1-PBS



Basic Information

Catalog Number: 60715-1-PBS	GenBank Accession Number: BC030155	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 4982	CloneNo.: 5F11E8
Source: Mouse	UNIPROT ID: O00300	
Isotype: IgG1	Full Name: tumor necrosis factor receptor superfamily, member 11b	
Immunogen Catalog Number: AG2139	Calculated MW: 401 aa, 46 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51028-1, TNFRSF11B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60715-1-PBS. Detection antibody: 60715-2-PBS. Standard:Ag2139. Range: 0.781-100 ng/mL