

For Research Use Only

# Transthyretin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60721-2-PBS



## Basic Information

Catalog Number:

60721-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

EG0931

GenBank Accession Number:

NM\_000371.4

GeneID (NCBI):

7276

UNIPROT ID:

P02766

Full Name:

transthyretin

Calculated MW:

16 kDa

Purification Method:

Protein G purification

CloneNo.:

3D4G2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51034-1, Transthyretin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60721-1-PBS. Detection antibody: 60721-2-PBS. Standard:Eg0931. Range: 0.098-100 ng/mL