

For Research Use Only

LGALS9/Galectin-9 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60731-1-PBS



Basic Information

Catalog Number:

60731-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC110340

GeneID (NCBI):

3965

UNIPROT ID:

O00182

Full Name:

lectin, galactoside-binding, soluble, 9

Calculated MW:

323 aa, 36 kDa

Purification Method:

Protein G purification

CloneNo.:

2G11B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51049-1, LGALS9/Galectin-9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60731-1-PBS. Detection antibody: 60731-2-PBS. Standard:Eg0084. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51049-2, LGALS9/Galectin-9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60731-1-PBS. Detection antibody: 60731-3-PBS. Standard:Eg0084. Range: 0.195-100 ng/mL