

For Research Use Only

CD83 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60744-4-PBS



Basic Information

Catalog Number:

60744-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27435

GenBank Accession Number:

BC030830

GeneID (NCBI):

9308

ENSEMBL Gene ID:

ENSG00000112149

UNIPROT ID:

Q01151

Full Name:

CD83 molecule

Calculated MW:

205 aa, 23 kDa

Purification Method:

Protein G purification

CloneNo.:

2F10B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51077-3, CD83 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60744-4-PBS. Detection antibody: 60744-2-PBS. Standard:Ag27435. Range: 0.195-100 ng/mL.