

For Research Use Only

Collagen Type XVIII Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60750-2-PBS



Basic Information

Catalog Number:

60750-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

EG0254

GenBank Accession Number:

NM_030582

GeneID (NCBI):

80781

UNIPROT ID:

P39060

Full Name:

collagen, type XVIII, alpha 1

Calculated MW:

178kd

Purification Method:

Protein A purification

CloneNo.:

1C8A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51199-1, COL18A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60750-1-PBS. Detection antibody: 60750-2-PBS. Standard:Eg0254. Range: 0.781-100 ng/mL