

For Research Use Only

FCN3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60753-1-PBS



Basic Information

Catalog Number:

60753-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30255

GenBank Accession Number:

BC020731

GeneID (NCBI):

8547

UNIPROT ID:

O75636

Full Name:

ficolin (collagen/fibrinogen domain containing) 3 (Hakata antigen)

Calculated MW:

33 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F7B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51200-2, Ficolin 3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60751-2-PBS. Detection antibody: 60753-1-PBS. Standard:Ag2449. Range: 6.25-200 ng/mL