

For Research Use Only

SURVIVIN-Specific Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60760-1-PBS



Basic Information

Catalog Number:

60760-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

GenBank Accession Number:

NM_001168

GeneID (NCBI):

332

UNIPROT ID:

O15392

Full Name:

baculoviral IAP repeat-containing 5

Calculated MW:

18 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

5G1B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51208-3, SURVIVIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66495-2-PBS. Detection antibody: 60760-1-PBS. Standard:Ag20958. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51208-4, SURVIVIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66495-4-PBS. Detection antibody: 60760-1-PBS. Standard:Ag20958. Range: 3.125-100 ng/mL