

For Research Use Only

# NOTCH1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60762-2-PBS



## Basic Information

Catalog Number:

60762-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG28356

GenBank Accession Number:

NM\_017617

GeneID (NCBI):

4851

UNIPROT ID:

P46531

Full Name:

Notch homolog 1, translocation-associated (Drosophila)

Calculated MW:

273 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G8B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51089-1, NOTCH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60762-1-PBS. Detection antibody: 60762-2-PBS. Standard:Ag28356. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51089-2, NOTCH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60762-3-PBS. Detection antibody: 60762-2-PBS. Standard:Ag28356. Range: 0.781-100 ng/mL.