

For Research Use Only

# CD180 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60768-5-PBS



## Basic Information

Catalog Number:

60768-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG28488

GenBank Accession Number:

BC109069

GeneID (NCBI):

4064

UNIPROT ID:

Q99467

Full Name:

CD180 molecule

Calculated MW:

661 aa, 74 kDa

Purification Method:

Protein A purification

CloneNo.:

2C6E2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51097-3, CD180 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60768-3-PBS. Detection antibody: 60768-5-PBS. Standard:Ag28488. Range: 0.391-50 ng/mL