

For Research Use Only

Cytokeratin 17 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60771-3-PBS



Basic Information

Catalog Number:

60771-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17657

GenBank Accession Number:

BC000159

GeneID (NCBI):

3872

UNIPROT ID:

Q04695

Full Name:

keratin 17

Calculated MW:

432 aa, 48 kDa

Purification Method:

Protein G purification

CloneNo.:

1H2A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51099-2, Cytokeratin 17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60771-3-PBS. Detection antibody: 60771-4-PBS. Standard:Ag17657. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP51099-4, Cytokeratin 17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60771-3-PBS. Detection antibody: 60771-6-PBS. Standard:Ag17657. Range: 1.563-200 ng/mL