

For Research Use Only

# MITF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60772-1-PBS



## Basic Information

Catalog Number: 60772-1-PBS	GenBank Accession Number: BC012503	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4286	CloneNo.: 1D1B5
Source: Mouse	UNIPROT ID: O75030	
Isotype: IgG2b	Full Name: microphthalmia-associated transcription factor	
Immunogen Catalog Number: AG14651	Calculated MW: 91 aa, 10 kDa, 59 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP51100-1, MITF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60772-1-PBS. Detection antibody: 60772-2-PBS. Standard:Ag14651. Range: 0.781-200 ng/mL.