

For Research Use Only

# LAMC3-Specific Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60773-1-PBS



## Basic Information

Catalog Number:

60773-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19466

GenBank Accession Number:

NM\_006059

GeneID (NCBI):

10319

UNIPROT ID:

Q9Y6N6

Full Name:

laminin, gamma 3

Calculated MW:

1575 aa, 171 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D7B12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51059-3, LAMC3-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60773-1-PBS. Detection antibody: 60773-2-PBS. Standard:Ag19466. Range: 0.781-12.5 ng/mL