

For Research Use Only

# SMARCA2-Specific Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60774-1-PBS



## Basic Information

Catalog Number:

60774-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33074

GenBank Accession Number:

NM\_003070

GeneID (NCBI):

6595

UNIPROT ID:

P51531

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

Calculated MW:

181 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F8E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51103-1, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-1-PBS. Detection antibody: 60774-2-PBS. Standard:Ag33074. Range: 3.125-100 ng/mL