

For Research Use Only

SMARCA2-Specific Monoclonal antibody, PBS Only (Detector)

Catalog Number:60774-2-PBS



Basic Information

Catalog Number:

60774-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33074

GenBank Accession Number:

NM_003070

GeneID (NCBI):

6595

UNIPROT ID:

P51531

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

Calculated MW:

181 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H1G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51103-1, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-1-PBS. Detection antibody: 60774-2-PBS. Standard:Ag33074. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51103-2, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-3-PBS. Detection antibody: 60774-2-PBS. Standard:Ag33074. Range: 0.391-100 ng/mL