

For Research Use Only

# VHL Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60776-4-PBS



## Basic Information

Catalog Number:

60776-4-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28860

GenBank Accession Number:

BC058831

GeneID (NCBI):

7428

UNIPROT ID:

P40337

Full Name:

von Hippel-Lindau tumor suppressor

Calculated MW:

172 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

3E3F3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51106-3, VHL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60776-4-PBS. Detection antibody: 60776-2-PBS. Standard: Ag28860. Range: 0.195-3.125 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 60776-4-Ig against Human VHL was performed. The affinity constant is below 1 pM.